


# Microscopy Journal club

Nikita Vladimirov, AG Preibisch

## **RAPID: Real-time image-based autofocus for all wide-field optical microscopy systems**


L. Silvestri, M. C. Müllenbroich, I. Costantini, A. P. Di Giovanna, L. Sacconi, and F. S. Pavone



# Autofocusing in microscopy (and photography)


## Common methods

### Contrast-based



- ✓ Image-based
- X Slow
- X User must define contrast function


### Triangulation-based



- X Based on reflective surface (~~SPIM~~)
- X Media R.I. should **not** match coverslip (for reflection)
- ✓ Fast
- ✓ Independent of the imaging


### Phase detecton...

## Contrast-based autofocus





<http://graphics.stanford.edu/courses/cs178-10/applets/autofocusCD.html>

## Contrast-based AF in SPIM system (Royer et al, *Nat. Biotech.*, 2016)




**Phase detection methods**






<http://graphics.stanford.edu/courses/cs178-10/applets/autofocusPD.html>

## Phase-detection autofocus in microscopy (some other methods)



**Parallax method:** Sun et al, *Nano Lett.*, 2009


- X Usable only in sparse samples (single-molecule conditions)
- ✓ Very accurate
- ✓ Allows sub-pixel localization
- X Slow



**Dual-cam Instantscope:**

Guo et al, *Biomed Optics Express*, 2015

- X Usable only in high light conditions
- ✓ Very accurate
- ✓ Very fast
- X Requires 2 extra CCDs




**Single-cam Instantscope**



Liao et al, *Biomed Optics Express*, 2016

- X Uses pinholes: only in high light conditions
- ✓ Very accurate
- ✓ Very fast
- ✓ Requires 1 extra CCD


Silvestri et al, 2017 setup:




- ✓ Accurate
- ✓ Fast
- ✓ Works in low-light conditions










## Cell culture




## Histology slice




## C.elegans tracking






SPIM, cleared mouse brain



Fixed focus




*RAPID*






SPIM, vasculature in mouse brain



SPIM, nuclear staining in mouse brain



## Disadvantages?

- Unlikely to work in low-light conditions (single-molecule) – using 10% of the collected light
- Not super-accurate, 0.7 of DOF

<i>Method</i>	<i>Accuracy (DOF)</i>	<i>Range (DOF)</i>	<i>Demonstrated time stability</i>	<i>Min S/N</i>	<i>Min S/B</i>	<i>Fluorescence</i>	<i>Dense samples</i>
<i>Parallax</i> <sup>17</sup>	~ 0.01	~ 8	minutes	~ 120	~1300	Yes	No
<i>Instantscope</i> <sup>18</sup>	~ 0.3	~ 300	N/A	N/A	N/A	No	Yes
Single-camera <i>Instantscope</i> <sup>19</sup>	~ 0.11	~ 25	N/A	N/A	N/A	No	Yes
<i>RAPID</i>	~ 0.7	~ 70	> 12 hours	~ 4.5	~ 1.7	Yes	Yes


Other options?

AF requires *information* from the sample

AF photons can be separated from data photons in various ways:

- in **time**: spend some imaging time (contrast-based methods, Parallax)
- in **photon budget**: sacrifice some data photons (Parallax, RAPID)
- in **spectrum**: separate the AF wavelength
  - IR beam reflected from glass (triangulation, e.g. ASI Crisp)
  - IR image of the sample (McGorthy et al, 2013)
  - time of flight for a laser pulse (LG smartphones)

Hybrid systems: split time and photon budget (modern Canon, Sony high-end cameras)



## Summary

The 'best' AF method depends on your application:

- Photon budget: how much light you can sacrifice?
- Speed: how fast you need to re-focus?
- Accuracy: what's your goal AF precision?
- Sample density: do you have points sources?
- Any reliable fiducials surfaces or beads?
- Can AF use a spare band of spectrum?